

IMPLANT DESCRIPTION	HEIGHT OF ABUTMENT TISSUE CUFF ABOVE THE IMPLANT							
	0 mm	.73 mm	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
<b>DENTSPLY FRIADENT (800.426.7836)</b>								
FRIALIT-2 & XIVE 3.4			8805	8806	8807	8808	8809	
FRIALIT-2 & XIVE 3.8			8810	8811	8812	8813	8814	
FRIALIT-2 & XIVE 4.5			8815	8816	8817	8818	8819	
FRIALIT 2- & XIVE 5.5			8139 (1.3mm)	8140	8141	8142	8143	
ANKYLOS				08064 (1.65mm)	8065	8066	8067	8068
<b>FMZ GmbH-ALPHATECH</b>								
Gingiformer Extra-Line Locator Implant Abutments (Available Exclusively from FMZ-GmbH Alphatech)								
<b>IMBIONIC IMPLANTAT</b>								
IMBIONIC Locator Implant Abutments (Available Exclusively from Imbionic Implantant)								
<b>IMPLANT DIRECT www.implantdirect.com</b>								
3.7 SCREW PLANT/LEGACY	8661		8662	8663	8664	8665	8625	8626
4.7 SCREW PLANT/LEGACY	8671		8672	8673	8674	8675	8627	8628
5.7 SCREW PLANT/LEGACY			8190	8191	8192	8193	8194	8195
3.5 REPLANT	8756		8757	8758	8759	8760	8734	8735
4.3 REPLANT	8761		8762	8763	8764	8765	8753	8754
5.0 REPLANT			8766	8767	8768	8769	8770	
<b>IMTEC CORPORATION (800.879.9799)</b>								
SELF TAPPING 3.75 (HEX)		8681		8682	8683	8684	8685	8696
UNIVERSAL 3.4, 4.0 (HEX)		8681		8682	8683	8684	8685	8696
MINI 3.3 (HEX)		8856 (.91 mm)		8857	8858	8859	8860	
ENDURE Locator Implant Abutments (Available Exclusively from Imtec Corporation)								
<b>INNOVA / SYBRON / ORALTRONICS (800.898.6261)</b>								
4.1 ENDOPORE (.7mm HEX)		8681		8682	8683	8684	8685	8696
3.25 ENTEGRA (.7mm HEX) 4.1 PLATFORM 03B		8681		8682	8683	8684	8685	8696
4.0 ENTEGRA (.7mm HEX)		8681		8682	8683	8684	8685	8696
3.5 ENDOPORE (.9mm HEX)			8851	8852	8853	8854	8855	
3.25 ENTEGRA (.9mm HEX) 3.5 PLATFORM 03M			8851	8852	8853	8854	8855	
4.1 ENDOPORE (.9mm HEX)			(There is no abutment available for this product)					
4.7 ENTEGRA (.9mm HEX)			8691	8692	8693	8694	8695	
5.0 ENDOPORE (.9mm HEX)			8691 (1.35mm)	8692	8693	8694	8695	
4.1, 5.0 ENDOPORE One-Stage (Internal Connection)		8950		8951	8952	8953	8954	8955
4.1 ENDOPORE (Internal Hex Connection)	8145		8146	8147	8148	8149	8150	
5.0 ENDOPORE (Internal Hex Connection)			8164	8165	8166	8167	8168	8169
3.25 / 3.3 Pitt-Easy			8311	8312	8313	8314	8315	8316
3.75 / 3.8 Pitt-Easy			8321	8322	8323	8324	8325	8326
4.0 / 4.1 Pitt-Easy			8331	8332	8333	8334	8335	8336
4.9 Pitt-Easy			8341	8342	8343	8344	8345	8346
<b>INTRA-LOCK SYSTEM INTERNATIONAL (877.866.0657)</b>								
3.75, 4.0 THREADED (HEX)		8681		8682	8683	8684	8685	8696
4.3 THREADED CONIC (HEX)			8781	8782	8783	8784	8785	
4.75 THREADED (HEX)			8727 (1.42 mm)	8728	8729	8730	8731	
5.0 THREADED CONIC (HEX)			8727 (1.42 mm)	8728	8729	8730	8731	
INTERNAL CONNECT Locator Implant Abutments (Available Exclusively from Intra-Lock Corporation)								
<b>KYOCERA</b>								
4.2 POI		8716 (.35mm)	8717	8718	8719	8720	8721	