

# LOCATOR F-Tx™ FIXED ATTACHMENT SYSTEM

## U.S. RETAIL PRODUCT ORDERING INFORMATION AND PRICE LIST\*

\*Products and pricing listed in this guide may not be available in all countries.

# LOCATOR F-Tx™ FIXED ATTACHMENT SYSTEM

Each LOCATOR F-Tx Abutment is delivered in convenient all-in-one packaging that includes all the sterile components needed for each implant.

**Order for your favorite implant systems using this product ordering and reference guide.**

LOCATOR F-Tx  
Abutment




**\$249 PRICE**  
INCLUDES:

- 1 - Abutment
- 1 - Denture Attachment Housing with Pre-inserted Processing Ball
- 1 - Low (**Blue**) Retention Ball
- 1 - Medium (**Tan**) Retention Ball
- 1 - High (**Green**) Retention Ball
- 1 - Additional Processing Ball
- 2 - Block-Out Spacer (Packaged separately)




## SELECTING THE PROPER CUFF HEIGHT FOR THE LOCATOR F-Tx ABUTMENT



1. Using a periodontal probe, measure from the top of the implant to the highest point of the gingiva.
2. The cuff size should match the tissue height, as this will place the prosthetic finish line right at the gingival height and maximizes the retention surface of the Denture Attachment Housing in the prosthesis. The cuff margin is at the widest portion of the sphere (midpoint of the Abutment sphere). For subgingival finish line, choose a shorter cuff size.

# LOCATOR F-Tx™ FIXED ATTACHMENT SYSTEM

## Abutment Ordering Information



**LOCATOR F-Tx ABUTMENTS ARE IMPLANT SPECIFIC**

**TO ORDER, YOU WILL NEED THE FOLLOWING INFORMATION:**

- 1.** Name of the Implant System
- 2.** Implant Diameter
- 3.** Connection type and platform diameter
- 4.** Cuff height selection criteria

IMPLANT COMPANY /DESCRIPTION	LOCATOR F-Tx ATTACHMENT SYSTEM (Includes Abutment)					
	CUFF HEIGHT					
	1.5mm	2mm	3mm	4mm	5mm	6mm

BIOHORIZONS						
<b>INTERNAL HEX CONNECTION</b> (Compatible with Tapered Plus, Tapered Internal Implant, Tapered Tissue Level, Single Stage & Internal Parallel Wall Implants)						
3.5mm Internal Hex Platform	10200-01	10200-02	10200-03	10200-04	10200-05	10200-06
4.5mm Internal Hex Platform	10201-01	10201-02	10201-03	10201-04	10201-05	10201-06

BIOMET 3i						
<b>CERTAIN INTERNAL CONNECTION</b> (Compatible with T3®, NanoTite® and OSSEOTITE® Implants)						
3.4mm Diameter CERTAIN Platform	10300-01	10300-02	10300-03	10300-04	10300-05	10300-06
4.1mm Diameter CERTAIN Platform	10301-01	10301-02	10301-03	10301-04	10301-05	10301-06
<b>EXTERNAL HEX CONNECTION</b> (Compatible with T3®, NanoTite® and OSSEOTITE® Implants)						
4.1mm Diameter External Hex Platform		10511-02	10511-03	10511-04	10511-05	10511-06

BLUE SKY BIO						
<b>BONE LEVEL CONNECTION</b> (Compatible with BIO Quattro Implants)						
Narrow Platform (NP)	10900-01	10900-02	10900-03	10900-04	10900-05	10900-06
Regular Platform (RP)	10901-01	10901-02	10901-03	10901-04	10901-05	10901-06
<b>TISSUE LEVEL CONNECTION</b> (Compatible with BIO One Stage Implants)						
Regular Platform (RP)		10906-02	10906-03	10906-04	10906-05	10906-06

DENTSPLY						
<b>ANKLYOS® CONICAL CONNECTION</b> (Compatible with Ankylos C/X Implants)						
Ankylos Implant Platform		10715-02	10715-03	10715-04	10715-05	10715-06
<b>ASTRA TECH EV CONNECTION</b> (Compatible with Astra Tech EV Implants)						
3.6mm Conical EV	10700-01	10700-02	10700-03	10700-04	10700-05	10700-06

# LOCATOR F-Tx™ FIXED ATTACHMENT SYSTEM

## Abutment Ordering Information (Continued)

IMPLANT COMPANY /DESCRIPTION	LOCATOR F-Tx ATTACHMENT SYSTEM (Includes Abutment)					
	CUFF HEIGHT					
	1.5mm	2mm	3mm	4mm	5mm	6mm
<b>ASTRA TECH EV CONNECTION</b> continued (Compatible with Astra Tech EV Implants)						
4.2mm Conical EV	10701-01	10701-02	10701-03	10701-04	10701-05	10701-06
<b>ASTRA TECH TX CONNECTION</b> (Compatible with Astra Tech TX Implants)						
ASTRA TX (ACQUA)	10706-01	10706-02	10706-03	10706-04	10706-05	10706-06
ASTRA TX (LILAC)	10707-01	10707-02	10707-03	10707-04	10707-05	10707-06
<b>NEODENT</b>						
<b>MORSE TAPER CONNECTION</b> (Compatible with CM Implants)						
CM Implant Platform		10715-02	10715-03	10715-04	10715-05	10715-06
<b>NOBEL BIOCARE</b>						
<b>CONICAL CONNECTION</b> (Compatible with NobelActive®, NobelParallel™ and NobelReplace® Implants)						
Narrow Platform (NP) Conical	10506-01	10506-02	10506-03	10506-04	10506-05	10506-06
Regular Platform (RP) Conical	10507-01	10507-02	10507-03	10507-04	10507-05	10507-06
<b>EXTERNAL HEX CONNECTION</b> (Compatible with Brånemark and NobelSpeedy™ Groovy Implants)						
Regular Platform (RP) External Hex		10511-02	10511-03	10511-04	10511-05	10511-06
<b>TRI-CHANNEL CONNECTION</b> (Compatible with NobelReplace, Replace Select™, NobelSpeedy™ Implants)						
Narrow Platform (NP) Tri-Channel	10500-01	10500-02	10500-03	10500-04	10500-05	10500-06
Regular Platform (RP) Tri-Channel	10501-01	10501-02	10501-03	10501-04	10501-05	10501-06
<b>SOUTHERN IMPLANTS</b>						
<b>TRI-NEX CONNECTION</b> (Compatible with Tri-Nex Co-Axis Implant, Tri-Nex Implant)						
3.5mm Tri-Nex Platform	10500-01	10500-02	10500-03	10500-04	10500-05	10500-06
4.3mm Tri-Nex Platform	10501-01	10501-02	10501-03	10501-04	10501-05	10501-06
<b>ZIMMER BIOMET</b>						
<b>INTERNAL HEX CONNECTION</b> (Compatible with Tapered Screw-Vent®, Screw-Vent®, Trabecular Metal™ Implants)						
3.5mm Diameter Platform	10200-01	10200-02	10200-03	10200-04	10200-05	10200-06
4.5mm Diameter Platform	10201-01	10201-02	10201-03	10201-04	10201-05	10201-06

# LOCATOR F-Tx™ FIXED ATTACHMENT SYSTEM

## Ancillary Components

PART #	DESCRIPTION		PRICE
<b>PROCESSING COMPONENTS</b>			
10041	Denture Attachment Housing (Qty. 2)		\$80.00
10032	Black Processing Ball (Qty. 2)		\$35.00
10033	Low Retention Ball ( <b>Blue</b> ), (Qty. 2)		\$35.00
10034	Medium Retention Ball ( <b>Tan</b> ), (Qty. 2)		\$35.00
10035	High Retention Ball ( <b>Green</b> ), (Qty. 2)		\$35.00
10036	Healing Cap (Qty. 2)		\$40.00
10037	Abutment Analog (Qty. 2)		\$40.00
10040	Waxing Cap (Qty. 2)		\$40.00
10043	Impression Coping (Qty. 2)		\$80.00
10030	Processing Cap (Qty. 2)		\$40.00
10039	Polishing Cap (Qty. 2)		\$40.00

# LOCATOR F-Tx™ FIXED ATTACHMENT SYSTEM

## Instrumentation and Prosthesis Removal Tools

PART #	DESCRIPTION	PRICE
<b>DRIVERS</b>		
10028	LOCATOR F-Tx Abutment Driver, Manual (used with square drive connections)	\$55.00
10029	LOCATOR F-Tx Abutment Driver, Latch (used with latch-lock connections)	\$55.00
10042	LOCATOR F-Tx Retention Ball Hex Driver	\$55.00
<b>PROSTHESIS REMOVAL — HYDRAULIC LIFT SYSTEM</b>		
10080	LOCATOR F-Tx Inflation Device	\$99.00
10045	LOCATOR F-Tx Prosthesis Booster Set, 10mm/8mm/8mm	\$49.99

10080



10045




2 x 8mm  
1 x 10mm


## PROSTHESIS REMOVAL UTILIZING THE HYDRAULIC LIFT SYSTEM



**1** Prosthesis Boosters inserted under the prosthesis between the prosthesis and soft tissue.



**2** Prosthesis Booster filled with water, and hydraulic displacement applied.



**3** Prosthesis disengages from the LOCATOR F-Tx Abutments.

# LOCATOR F-Tx™ FIXED ATTACHMENT SYSTEM

## Instrumentation and Prosthesis Removal Tools (Continued)

PART #	DESCRIPTION	PRICE
PROSTHESIS REMOVAL TOOL BAR AND LOOP		
10024	LOCATOR F-Tx Removal Tool Bar	\$99.00
10023	LOCATOR F-Tx Removal Loop	\$30.00




10023




10024


### PROSTHESIS REMOVAL UTILIZING THE TOOL BAR AND LOOP



**1** Place the Loop under the posterior portion of one side of the prosthesis.



**2** Slide the Bar handle through the two loops and place the Delrin tip on top of the contralateral posterior occlusal surface.



**3** **Gently** load the bar handle until one side of the prosthesis disengages from the LOCATOR F-Tx Abutments.



**LOCATOR F-Tx™**  
FIXED ATTACHMENT SYSTEM

VISIT

**ZESTDENT.COM**

24 HOURS A DAY | 7 DAYS A WEEK

TO PLACE ORDERS ONLINE OR TO FIND THE LATEST  
SUPPORT TOOLS FOR LOCATOR F-Tx™



ZEST | DANVILLE MATERIALS | PERIOSCOPY

ZEST DENTAL SOLUTIONS  
2875 LOKER AVENUE EAST  
CARLSBAD, CA 92010  
USA

TEL: 800.262.2310  
FAX: 760.743.7975  
EMAIL: SALES@ZESTDENT.COM  
**WWW.ZESTDENT.COM**

DISTRIBUTED BY: